

Adult Human Ex-Vivo Models for Preclinical Cardiac Safety Assessment of Drugs




Dr. Najah Abi-Gerges

VP of R&D

najah.abigergeres@anabios.com

Current Preclinical Cardiac Safety Approaches

Have Significant Limitations






Human cardiac models are urgently needed

for detection of drug-induced cardiotoxicity

Enabling Drug Discovery in Human Healthy and Diseased Tissues




- Tissue harvesting methods and solutions are designed to avoid ischemic damage and reperfusion injury
- Complete chain of custody, processing methods and rigorous QC ensure excellent tissue quality
- Large U.S.A.-based network ensures the availability of samples
- Excellent heart quality permits integrated human cardiac drug discovery at the preclinical stages

Human Ex-Vivo

Cardiac Safety Assessment at AnaBios

CELL-BASED ASSAYS

(Optimization of drugs)



>1000 **ex vivo**
human hearts tested
2-4 hearts / week

Arrhythmia & Inotropy
Ventricular or atrial myocytes contractility

Intracellular Ca²⁺ Dynamics
Ventricular Myocytes

Ion Channel Block
Ventricular or atrial myocytes V-clamp

Action Potential
Ventricular & Atrial Myocytes I-clamp

Cardiac Fibrosis
Cardiac Fibroblasts

TISSUE-BASED ASSAYS

(Nomination of drugs)

Pro-arrhythmia
*Action Potential
Ventricular Trabeculae*

Inotropy
*Contractility
Ventricular & Atrial Trabeculae*



Chronotropy
Spontaneous Action Potential Sinoatrial Node

Vaso-constriction Dilation
Coronary Rings

New Isolation Method Provides High Yield of Cardiomyocytes




1000 **ex vivo**
human hearts
tested
2-4 hearts / week



Non-Invasive Measurement of Contraction


Full Retention of Cardiomyocyte Functionality



- Bright-field imaging
- Low technical complexity
- No cytotoxic fluorescent reagents
- High information content


IonOptix: Sarcomere shortening measured by digital cell geometry tracking; stimulation frequency 1Hz

Strong Correlation Between Electrical and Mechanical Abnormalities in the Human Heart




Nador *et al.*, (1991) *Circulation* 84:1530-1542 and other papers (De Ferrari *et al.*, 1994; Nakayama *et al.*, 1998; Haugaa *et al.*, 2009; Ferrari & Schwartz, 2009; Belardinelli *et al.*, 2009)

Markers of Pro-arrhythmia and Contractility Risk



Normal condition





Drug treatment

- After-contraction represents the mechanical manifestation of triggered EAD
- TdP arrhythmia arises from PVCs due to triggered EADs (Kaumann et al. 1968; Noda et al. 2014)

Validating Clinical Relevance


- Validated 33 clinical well characterized controls (24 CiPA / 32 JiCSA):
 - 1) 23 pro-arrhythmic drugs
 - 2) 10 non-pro-arrhythmic drugs
 - 3) Each drug was tested at multiples of the free Effective Therapeutic Plasma Concentration (fETPC, mimic pharmacokinetic aspect)
 - 4) Establish pharmacodynamic exposure response

AC for Predicting Pro-arrhythmic Risk



Prediction of Pro-arrhythmic Risk of Drugs with 96% sensitivity

Table 1. Pro-arrhythmia prediction of the adult human primary cardiomyocyte-based model



| Drug name | Clinical TdP risk | Pro-arrhythmia risk at 10-fold FETPC | | | | | |
|-----------------------------|-------------------|---|--|--|--|---|---|
| | | ANABIOS Adult human primary ventricular (Sarc. Short., AC) Nguyen et al., 2017 | AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA FPD) Qu et al., 2015 | AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA EAD) Qu et al., 2016 | JiCSA hiPSC-derived cardiomyocytes (iCell®, MEA Score) Ando et al., 2017 | FDA hiPSC-derived cardiomyocytes (iCell®, MEA Arrhythmia) Blinova et al., 2017 | FDA hiPSC-derived cardiomyocytes (Cor.4U, MEA Arrhythmia) Blinova et al., 2017 |
| Ajmaline | | | Not tested | Not tested | | Not tested | Not tested |
| Astemizole ^a | | False negative | Not tested | Not tested | | Not tested | Not tested |
| Azimilide ^a | | | Not tested | Not tested | Not tested | Not tested | Not tested |
| Bepiridil ^a | | | Not tested | Not tested | False negative | False negative | False negative |
| Chlorpromazine ^a | | | Not tested | Not tested | False negative | False negative | False negative |
| Cisapride ^a | | | | False negative | | False negative | False negative |
| Clarithromycin ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Clozapine ^a | | | Not tested | Not tested | False negative | Not tested | Not tested |
| D, L-Sotalol ^a | | | | | | Not tested | Not tested |
| Disopyramide ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Dofetilide ^a | | | | | | | |
| Domperidone ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Droperidol ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Erythromycin | | | Not tested | Not tested | | Not tested | Not tested |
| Flecainide | | | | | | Not tested | Not tested |
| Ibutilide ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Moxifloxacin | | | | Not tested | | False negative | False negative |
| Ondansetron ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Procainamide | | | Not tested | Not tested | | Not tested | Not tested |
| Quinidine ^a | | | Not tested | Not tested | | | |
| Sematilide | | | Not tested | Not tested | | Not tested | Not tested |
| Terodiline | | | False negative | False negative | | Not tested | Not tested |
| Vandetanib ^a | | | Not tested | Not tested | | Not tested | Not tested |

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: negative pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell® hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; FPD: Field Potential Duration; JiCSA: Japan iPS Cardiac Safety Assessment; FDA: Food and Drug Administration; Cor.4U: hiPSC-derived cardiomyocytes from AxioGenesis AG; EAD: Early afterdepolarization.

Prediction of Non-Pro-arrhythmic Drugs with 100% Specificity




Table 3. Pro-arrhythmia prediction of the adult human primary cardiomyocyte-based model

| Drug name | Clinical TdP risk | Pro-arrhythmia risk at 10-fold fETPC | | | | | |
|---------------------------|-------------------|--|--|---|---|---|---|
| | | ANABIOS primary ventricular (Sarc. short AC) Nguyen et al., 2017 | AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA FPD) Qu et al., 2015 | AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA EAD) Qu et al., 2016 | JiCSA hiPSC-derived cardiomyocytes (iCell®, MEA Score) Ando et al., 2017 | FDA hiPSC-derived cardiomyocytes (iCell®, MEA Arrhythmia) Blinova et al., 2017 | FDA hiPSC-derived cardiomyocytes (Cor.4U, MEA Arrhythmia) Blinova et al., 2017 |
| Diltiazem ^a | | | Not tested | Not tested | | | |
| Diphenhydramine | | | Not tested | Not tested | False positive | Not tested | Not tested |
| Loratadine ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Mexiletine ^a | | | False positive | Not tested | False positive | Quiescent | |
| Mibefradil | | | Not tested | Not tested | | | |
| Nifedipine ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Nitrendipine ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Ranolazine ^a | | | False positive | | False positive | False positive | |
| Tamoxifen ^a | | | Not tested | Not tested | | Not tested | Not tested |
| Verapamil ^a | | | Not tested | Not tested | | | Quiescent |

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: negative pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell® hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; FPD: Field Potential Duration; JiCSA: Japan iPS Cardiac Safety Assessment; FDA: Food and Drug Administration; Cor.4U: hiPSC-derived cardiomyocytes from Axiogenesis AG; EAD: Early afterdepolarization.

ICH S7B IWG Recognizes the Value of Human Primary Cardiomyocytes for Pro-arrhythmia Assessment



Final Concept Paper

ICH S7B and E14 Q&A


Endorsed by the ICH Assembly on 15 November 2018

ICH S7B recommends Follow-up Studies (Sec. 2.3.5) to inform the integrated risk assessment if a test article blocks the hK_v11.1 IKr current (hERG) or prolongs the QT interval. These could include the test article effects on additional ionic currents, and the use of *in vitro* and *in vivo* assays. Newer assays and technologies such as *in silico* ventricular models, and **human primary** and induced pluripotent stem cell-derived cardiomyocytes, can provide insights into the relative proarrhythmic liability of test articles. Guidance is needed on when and how these novel approaches play a role in determining the proarrhythmic risk to inform clinical development.

Validating Clinical Relevance of Negative Inotropes

- Validated 33 clinical well characterized controls:
 - 1) 27 multichannel blockers (mainly K^+ , Na^+ and Ca^{2+} channels) as positive controls
 - 2) 6 selective hERG blockers as negative controls

Verapamil Induces Negative Inotropic Effect



Identification of Negative Inotropic Effects and Determination of Exposure Responses





Table 5. Sarcomere shortening effects for reference drugs measured in adult human primary cardiomyocytes


| Drug name | Top test concentration (μM) | Human myocyte effect | IC ₅₀ (μM) | Ratio (IC ₅₀ /fETPC) |
|-----------------------------|-----------------------------|----------------------|-----------------------|---------------------------------|
| Ajmaline | 1.95 | -ve inotrope | 2 | 31 |
| Astemizole ^a | 0.009 | No effect | >0.009 | 30 |
| Azimilide ^a | 2.1 | -ve inotrope | 1.07 | 15 |
| Bepridil ^a | 0.96 | -ve inotrope | 0.7 | 22 |
| Chlorpromazine ^a | 1.04 | -ve inotrope | 1.02 | 28 |
| Cisapride ^a | 0.26 | -ve inotrope | 0.02 | 8 |
| Clarithromycin ^a | 120 | -ve inotrope | 16 | 13 |
| Clozapine ^a | 2.13 | -ve inotrope | 1.5 | 21 |
| D, L-Sotalol ^a | 450 | No effect | >450 | >30 |
| Disopyramide ^a | 21 | -ve inotrope | 9.3 | 13 |
| Dofetilide ^a | 0.2 | No effect | >0.2 | >100 |
| Domperidone ^a | 2 | -ve inotrope | 0.2 | 10 |
| Droperidol ^a | 0.48 | -ve inotrope | 0.18 | 11 |
| Erythromycin | 5.1 | No effect | >5.1 | >30 |
| Flecainide | 22.6 | -ve inotrope | 1.1 | 2 |
| Ibutilide ^a | 3 | -ve inotrope | 2 | 20 |
| Moxifloxacin | 329 | No effect | >329 | >30 |
| Ondansetron ^a | 11.2 | -ve inotrope | 14 | 34 |
| Procainamide | 1625 | -ve inotrope | 2215 | 38 |
| Quinidine ^a | 100 | -ve inotrope | 3.6 | 1 |
| Sematilide | 133 | No effect | >133 | >30 |
| Terodiline | 4.35 | -ve inotrope | 0.7 | 5 |
| Vandetanib ^a | 9 | -ve inotrope | 2.7 | 9 |
| Diltiazem ^a | 3.84 | -ve inotrope | 1 | 8 |
| Diphenhydramine | 1.02 | -ve inotrope | 0.6 | 17 |
| Loratadine ^a | 0.0135 | -ve inotrope | 0.0175 | 35 |
| Mexiletine ^a | 75 | -ve inotrope | 0.9 | 0.4 |
| Mibefradil | 0.36 | -ve inotrope | 0.18 | 13 |
| Nifedipine ^a | 0.23 | -ve inotrope | 0.04 | 5 |
| Nitrendipine ^a | 0.091 | -ve inotrope | 0.06 | 18 |
| Ranolazine ^a | 200 | -ve inotrope | 17 | 9 |
| Tamoxifen ^a | 0.663 | -ve inotrope | 0.99 | 36 |
| Verapamil ^a | 10 | -ve inotrope | 0.04 | 2 |

IC₅₀: Concentration inducing 50% decrease in sarcomere shortening; Hill equation using SigmaPlot v13 was fitted to sarcomere shortening concentration-effect curves, assuming drugs would eventually cause complete inhibition of the contractility when they decreased sarcomere shortening by ≥25%. ^a: CiPA-selected drug; fETPC, free effective therapeutic plasma concentration.

Nguyen et al., 2017 FiP

Inhibition of Kinase Activity to Control Tumor Growth


Can Lead to Cardiotoxicity




Validating Clinical Relevance of Cancer Agents

- Validated 9 clinical well characterized controls :
 - 1) 4 toxic TKIs (Sorafenib, Vandetanib, AZD7762, Imatinib)
 - 2) 4 non-toxic TKIs (Erlotinib, Dasatinib, Afatinib, Gefitinib)
 - 3) One toxic anthracycline (Doxorubicin)
 - 4) Each drug was tested at multiples of the Cmax
 - 5) Each concentration was perfused for 5 mins

Sorafenib Induces Functional Cardiotoxicity




Afatinib Induces No Functional or Structural Cardiotoxicity



$0.3\mu M = 3\text{-fold } C_{max}$
 $1\mu M = 10\text{-fold } C_{max}$
 $3\mu M = 30\text{-fold } C_{max}$

Tyrosine Kinase Inhibitors Affect Human Cardiomyocyte Contractility







| TKI | Clinical contractility risk | Human cardiomyocyte contractility | C _{max} (μM) | IC ₅₀ (μM) | Ratio (IC ₅₀ /C _{max}) |
|------------|-----------------------------|-----------------------------------|-----------------------|-----------------------|---|
| Sorafenib | | | 3.4 | 1.2 | 0.35 |
| Vandetanib | | | 1.8 | 4.6 | 2.55 |
| AZD7762 | | | 0.12 | 0.8 | 6 |
| Imatinib | | | 5 | 44 | 8 |
| Erlotinib | | | 2.5 | >10* | >4 |
| Dasatinib | | | 0.01 | >0.3 | >30 |
| Afatinib | | | 0.1 | >3 | >30 |
| Gefitinib | | | 0.1 | >3 | >30 |


Similar human cardiac tissue data recently published by Schneider C et al., 2018 Nature Scientific Reports

*: Limit of solubility

Low Inter- and Intra-Heart Variability






Human cardiomyocyte Contractility Model Differentiates Pro-arrhythmic from Non-Pro-arrhythmic TKIs



| TKI | Clinical contractility risk | Human cardiomyocyte pro-A | Cmax (μM) | Pro-A conc. (μM) | Ratio (Pro-A conc./C _{max}) |
|------------|-----------------------------|---------------------------|------------------------|-------------------------------|---------------------------------------|
| Vandetanib | | | 1.8 | 3 | 1.6 |
| AZD7762 | | | 0.12 | 0.3 | 2.5 |
| Imatinib | | | 5 | 6 | 1.2 |
| Sorafenib | | | 3.4 | >10 | >3 |
| Erlotinib | | | 2.5 | >10* | >4 |
| Dasatinib | | | 0.01 | >0.3 | >30 |
| Afatinib | | | 0.1 | >3 | >30 |
| Gefitinib | | | 0.1 | >3 | >30 |

*: Limit of solubility

Doxorubicin, Anthracycline Agent, Affects Human Cardiomyocyte Contractility and Induces Pro-arrhythmia



Adult Human Heart Kinome Profiling

**Cancer-driving
kinases**



Expression and role
of kinases in cancer
are well understood

Heart kinases




- Expression & function of kinases in cardiac tissue are poorly characterized
- Mechanisms of KI-induced cardiotoxicity are not fully understood

Positive Inotropy Assessment - Validation Set

Targets 12 Different Mechanisms of Action

| Mechanism of Action | Drug |
|--|--------------------|
| Na ⁺ /K ⁺ pump inhibition | Digoxin |
| Na ⁺ /K ⁺ pump inhibition | Ouabain |
| Na ⁺ /Ca ²⁺ exchanger inhibition | SEA-0400 |
| Myosin activation | Omecamtiv Mecarbil |
| Myosin activation | EMD-57003 |
| Ca ²⁺ sensitization | Levosimendan |
| Non-selective b-adrenoceptor activation | Isoproterenol |
| Non-selective b-adrenoceptor activation | Epinephrine |
| b1-adrenoceptor activation | Dobutamine |
| PDE3 inhibition | Milrinone |
| PDE inhibition | IBMX |
| Ca ²⁺ channel activation | Bay-K 8644 |
| Adenylyl cyclase activation | Forskolin |
| Adenylyl cyclase activation | NKH-477 |
| SERCA activation | N106 |
| RyR activation | Caffeine |










Isoproterenol Induces Positive Inotropic Effect




Baseline level

Iso level


Identification of Positive Inotropic Effects and Determination of Exposure Responses



Differential Effects of Positive Inotropes







Cluster Analysis is Used to Mechanistically Fingerprint Compounds with Inotropic Effects



Heatmap data generated from 4th concentration data. Red and green colors indicate decrease and increase of >25% and 10% change, respectively. Black colors indicate no effect (<-25% < % change < 10%). Numbers in boxes indicate means % change relative to vehicle.


Segregation of Ca^{2+} -Dependent Mechanisms



2D PCA generated from top test concentration data. Blue and red colors indicate increase and decrease in Ca^{2+} , respectively. Ellipses show confidence intervals of 0.75.



Human Adult Cardiomyocyte Model

Integrated Into the CiPA Paradigm



This new potential approach was presented to the ICH S7B & E14 Discussion Group Meeting (Charlotte, NC, November 13, 2018).

MyoBLAZER™: Proprietary Technology Used to Measure the Biomarker Currently Undergoing Qualification with FDA



Isolation of human cardiac cells





Proprietary cell contractility testing station

Proprietary acquisition and analysis algorithms allow to establish drug effects with high level of predictivity of the human clinical response



Adult Human Cardiomyocyte Model

Early Primary Screening Tool



Adult Human Cardiomyocyte Model

Early Primary Cardiotoxicity Screening Tool

- Permit rigorous and integrated human cardiac drug discovery at the preclinical stages
- Differentiate safe from cardiotoxic drugs
- Can enable mechanistic assessment
- Predictive of clinical outcomes

Thank You!