

Higher Throughput Adult Human Primary Cardiomyocyte Model for Drug Safety Screening




Dr. Najah Abi-Gerges

VP of R&D

najah.abigergeres@anabios.com

Current Preclinical Cardiac Safety Approaches

Have Significant Limitations






Human cardiac models are urgently needed

for detection of drug-induced cardiotoxicity

Enabling Drug Discovery in Human Healthy and Diseased Tissues



- Tissue harvesting methods and solutions are designed to avoid ischemic damage and reperfusion injury
- Complete chain of custody, processing methods and rigorous QC ensure excellent tissue quality
- Large U.S.A.-based network ensures the availability of samples
- Excellent heart quality permits integrated human cardiac drug discovery at the preclinical stages

Human Ex-Vivo


Cardiac Safety Assessment at AnaBios

CELL-BASED ASSAYS

(Optimization of drugs)

TISSUE-BASED ASSAYS

(Nomination of drugs)



>1400 **ex vivo**
human hearts tested
2-4 hearts / week

Arrhythmia & Inotropy
Ventricular or atrial myocytes contractility

Intracellular Ca²⁺ Dynamics
Atrial or Ventricular Myocytes

Ion Channel Block
Ventricular or atrial myocytes V-clamp

Action Potential
Ventricular & Atrial Myocytes I-clamp

Cardiac Fibrosis
Cardiac Fibroblasts

Pro-arrhythmia
*Action Potential
Ventricular Trabeculae*

Inotropy
*Contractility
Ventricular & Atrial Trabeculae*



Chronotropy
Spontaneous Action Potential Sinoatrial Node

Vaso-constriction Dilation
Coronary Rings

New Isolation Method Provides High Yield of Adult Human Primary Cardiomyocytes




>1400 **ex vivo**
human hearts tested
2-4 hearts / week



Non-Invasive Measurement of Contraction


Full Retention of Cardiomyocyte Functionality



- Bright-field imaging
- Low technical complexity
- No cytotoxic fluorescent reagents
- High information content


IonOptix: Sarcomere shortening measured by digital cell geometry tracking; stimulation frequency 1Hz

Strong Correlation Between Electrical and Mechanical Abnormalities in the Human Heart




Nador *et al.*, (1991) *Circulation* 84:1530-1542 and other papers (De Ferrari *et al.*, 1994; Nakayama *et al.*, 1998; Haugaa *et al.*, 2009; Ferrari & Schwartz, 2009; Belardinelli *et al.*, 2009)

Strong Correlation Between Electrical and Mechanical Functions in the Adult Human Primary Cardiomyocyte




Baskin et al., 2016 JCI

Markers of Pro-arrhythmia and Contractility Risk



Normal condition




Drug treatment

- After-contraction represents the mechanical manifestation of triggered EAD
- TdP arrhythmia arises from PVCs due to triggered EADs (Kaumann et al. 1968; Noda et al. 2014)

Validating Clinical Relevance


- Validated 33 clinical well characterized controls (24 CiPA / 32 JiCSA):
 - 1) 23 pro-arrhythmic drugs
 - 2) 10 non-pro-arrhythmic drugs
 - 3) Each drug was tested at multiples of the free Effective Therapeutic Plasma Concentration (fETPC, mimic pharmacokinetic aspect)
 - 4) Establish pharmacodynamic exposure response

AC for Predicting Pro-arrhythmic Risk



Prediction of Pro-arrhythmic Risk of Drugs with 96% sensitivity

Table 1. Pro-arrhythmia prediction of the adult human primary cardiomyocyte-based model



Drug name	Clinical TdP risk	Pro-arrhythmia risk at 10-fold FETPC					
		ANABIOS Adult human primary ventricular (Sarc. Short., AC) Nguyen et al., 2017	AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA FPD) Qu et al., 2015	AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA EAD) Qu et al., 2016	JiCSA hiPSC-derived cardiomyocytes (iCell®, MEA Score) Ando et al., 2017	FDA hiPSC-derived cardiomyocytes (iCell®, MEA Arrhythmia) Blinova et al., 2017	FDA hiPSC-derived cardiomyocytes (Cor.4U, MEA Arrhythmia) Blinova et al., 2017
Ajmaline			Not tested	Not tested		Not tested	Not tested
Astemizole ^a		False negative	Not tested	Not tested		Not tested	Not tested
Azimilide ^a			Not tested	Not tested	Not tested	Not tested	Not tested
Bepiridil ^a			Not tested	Not tested	False negative	False negative	False negative
Chlorpromazine ^a			Not tested	Not tested	False negative	False negative	False negative
Cisapride ^a				False negative		False negative	False negative
Clarithromycin ^a			Not tested	Not tested		Not tested	Not tested
Clozapine ^a			Not tested	Not tested	False negative	Not tested	Not tested
D, L-Sotalol ^a						Not tested	Not tested
Disopyramide ^a			Not tested	Not tested		Not tested	Not tested
Dofetilide ^a							
Domperidone ^a			Not tested	Not tested		Not tested	Not tested
Droperidol ^a			Not tested	Not tested		Not tested	Not tested
Erthromycin			Not tested	Not tested		Not tested	Not tested
Flecainide						Not tested	Not tested
Ibutilide ^a			Not tested	Not tested		Not tested	Not tested
Moxifloxacin				Not tested		False negative	False negative
Ondansetron ^a			Not tested	Not tested		Not tested	Not tested
Procainamide			Not tested	Not tested		Not tested	Not tested
Quinidine ^a			Not tested	Not tested			
Sematilide			Not tested	Not tested		Not tested	Not tested
Terodiline			False negative	False negative		Not tested	Not tested
Vandetanib ^a			Not tested	Not tested		Not tested	Not tested

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: negative pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell® hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; FPD: Field Potential Duration; JiCSA: Japan iPS Cardiac Safety Assessment; FDA: Food and Drug Administration; Cor.4U: hiPSC-derived cardiomyocytes from AxioGenesis AG; EAD: Early afterdepolarization.

Human Stem Cell Cardiomyocytes Lack Phenotypic Stability for the Prediction of Pro-arrhythmia Risk

Table 2. Pro-arrhythmia prediction: Adult human primary cardiomyocytes *versus* stem cell-derived cardiomyocytes

Drug name	Clinical TdP risk	Pro-arrhythmia risk at 10-fold fETPC															
		ANABIOS Adult human primary ventricular (Sarc. Short. AC)	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 1, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 2, CLY AP, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 3, MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 4, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 5, MCS MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 1, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 2, CLY AP, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 3, MCS MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 4, ECR MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 5, MCS MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 6, ECR MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 7, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 8, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 9, AXN MEA, Arrhythmia) Blinova et al..	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 10, AMD MEA, Arrhythmia) Blinova et al..
Astemizole ^a		False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
Azimilide ^a							False negative										
Bepidil ^a			False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
Chlorpromazine ^a			False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
Cisapride ^a			False negative	False negative		False negative	False negative		False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
Clarithromycin ^a			False negative	False negative	False negative	False negative	False negative	False negative		False negative	False negative	False negative	False negative				False negative
Clozapine ^a			False negative	False negative	Quiescence	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
D, L-Sotalol ^a					False negative												
Disopyramide ^a			False negative	False negative		False negative	False negative			False negative	False negative	False negative				False negative	
Dofetilide ^a						False negative	False negative					False negative					
Domperidone ^a							False negative					False negative	False negative				
Droperidol ^a			False negative	False negative	False negative	False negative	False negative		False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative	False negative
Ibutilide ^a																	
Ondansetron ^a			False negative	False negative		False negative	False negative			False negative		False negative	False negative			False negative	
Quinidine ^a										False negative							
Vandetanib ^a			False negative	False negative	False negative	False negative	False negative				False negative				False negative		False negative

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: False negative non-pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell2 hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; AXN: Axion Biosystems; CLY: Clyde Biosciences; MCS: Multichannel Systems; ECR: ACEA CardioECR; AMD: Alpha MED Scientific.

Prediction of Non-Pro-arrhythmic Drugs with 100% Specificity




Table 3. Pro-arrhythmia prediction of the adult human primary cardiomyocyte-based model

Drug name	Clinical TdP risk	Pro-arrhythmia risk at 10-fold fETPC					
		ANABIOS primary ventricular (Sarc. short AC) Nguyen et al., 2017	AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA FPD) Qu et al., 2015	AMGEN hiPSC-derived cardiomyocytes (iCell®, MEA EAD) Qu et al., 2016	JiCSA hiPSC-derived cardiomyocytes (iCell®, MEA Score) Ando et al., 2017	FDA hiPSC-derived cardiomyocytes (iCell®, MEA Arrhythmia) Blinova et al., 2017	FDA hiPSC-derived cardiomyocytes (Cor.4U, MEA Arrhythmia) Blinova et al., 2017
Diltiazem ^a			Not tested	Not tested			
Diphenhydramine			Not tested	Not tested	False positive	Not tested	Not tested
Loratadine ^a			Not tested	Not tested		Not tested	Not tested
Mexiletine ^a			False positive	Not tested	False positive	Quiescent	
Mibefradil			Not tested	Not tested			
Nifedipine ^a			Not tested	Not tested		Not tested	Not tested
Nitrendipine ^a			Not tested	Not tested		Not tested	Not tested
Ranolazine ^a			False positive		False positive	False positive	
Tamoxifen ^a			Not tested	Not tested		Not tested	Not tested
Verapamil ^a			Not tested	Not tested			Quiescent

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: negative pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell® hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; FPD: Field Potential Duration; JiCSA: Japan iPS Cardiac Safety Assessment; FDA: Food and Drug Administration; Cor.4U: hiPSC-derived cardiomyocytes from Axiogenesis AG; EAD: Early afterdepolarization.

Human Stem Cell Cardiomyocytes Lack Phenotypic Stability for the Prediction of Non-Pro-arrhythmic Drugs

Table 4. Pro-arrhythmia prediction: Adult human primary cardiomyocytes *versus* stem cell-derived cardiomyocytes

Drug name	Clinical TdP risk	Pro-arrhythmia risk at 10-fold fETPC															
		ANABIOS Adult human primary ventricular cardiomyocytes (Sarc short AC) Nguyen et al., 2017	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 1, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 2, CLY AP, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 3, MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 4, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (Ncardia, Site 5, MCS MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 1, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 2, CLY AP, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 3, MCS MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 4, ECR MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 5, MCS MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 6, ECR MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 7, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 8, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 9, AXN MEA, Arrhythmia) Blinova et al., 2018	CiPA hiPSC-derived cardiomyocytes (iCell2, Site 10, AMD MEA, Arrhythmia) Blinova et al., 2018
Diltiazem ^a				Quiescence	Quiescence												
Loratadine ^a																	
Mexiletine ^a				Quiescence			Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence	Quiescence
Nifedipine ^a				Quiescence													
Nitrendipine ^a					Quiescence												
Ranolazine ^a				Quiescence	Quiescence		Quiescence				Quiescence	Quiescence	Quiescence	Quiescence	Quiescence		Quiescence
Tamoxifen ^a																	
Verapamil ^a				Quiescence	Quiescence	Quiescence					Quiescence		Quiescence	Quiescence		Quiescence	

^a: CiPA-selected drug; Red: positive pro-arrhythmia risk; Green: False negative non-pro-arrhythmia risk; hiPSC: human induced pluripotent stem cell (hiPSC); iCell2 hiPSC-derived cardiomyocytes from Cellular Dynamics; MEA: micro-electrode array; AXN: Axion Biosystems; CLY: Clyde Biosciences;

- Inability to demonstrate exposure response (High incidence of Quiescence)

ICH S7B IWG Recognizes the Value of Adult Human Primary Cardiomyocytes for Pro-arrhythmia Assessment



Final Concept Paper


ICH S7B and E14 Q&A

Endorsed by the ICH Assembly on 15 November 2018

ICH S7B recommends Follow-up Studies (Sec. 2.3.5) to inform the integrated risk assessment if a test article blocks the hK_v11.1 IKr current (hERG) or prolongs the QT interval. These could include the test article effects on additional ionic currents, and the use of *in vitro* and *in vivo* assays. Newer assays and technologies such as *in silico* ventricular models, and **human primary** and induced pluripotent stem cell-derived cardiomyocytes, can provide insights into the relative proarrhythmic liability of test articles. Guidance is needed on when and how these novel approaches play a role in determining the proarrhythmic risk to inform clinical development.

FDA Awards AnaBios Grant

[Resources](#) [Blog](#) [Journalists](#) [Log In](#) [Sign Up](#) [GDPR](#) [Send a Release](#)

 [News](#) [Products](#) [Contact](#)

[News in Focus](#) [Business & Money](#) [Science & Tech](#) [Lifestyle & Health](#) [Policy & Public Interest](#) [People & Culture](#)

FDA Awards AnaBios Grant to Further Develop Preclinical Assay Using Human Primary Cardiomyocytes

NEWS PROVIDED BY
[AnaBios](#) →
Sep 17, 2019, 08:00 ET


SHARE THIS ARTICLE
[f](#) [t](#) [in](#) [p](#) [e](#) [p](#)

SAN DIEGO, Sept. 17, 2019 /PRNewswire/ -- The U.S. Food and Drug Administration (FDA) has named AnaBios the recipient of a grant to further develop its unique drug discovery platform utilizing muscle cells from the human heart (cardiomyocytes). Funds from the grant will be used to develop a preclinical biomarker to identify the pro-arrhythmia risk of potential drugs based on contractility measurements in human adult primary cardiomyocytes. AnaBios' cardiomyocyte research has yielded highly translational information about cellular properties and drug-induced variations in cardiac function.

Validating Clinical Relevance of Negative Inotropes

- Validated 33 clinical well characterized controls:
 - 1) 27 multichannel blockers (mainly K^+ , Na^+ and Ca^{2+} channels) as positive controls
 - 2) 6 selective hERG blockers as negative controls

Verapamil Induces Negative Inotropic Effect



Identification of Negative Inotropic Effects and Determination of Exposure Responses





Table 5. Sarcomere shortening effects for reference drugs measured in adult human primary cardiomyocytes


Drug name	Top test concentration (μM)	Human myocyte effect	IC ₅₀ (μM)	Ratio (IC ₅₀ /fETPC)
Ajmaline	1.95	-ve inotrope	2	31
Astemizole ^a	0.009	No effect	>0.009	30
Azimilide ^a	2.1	-ve inotrope	1.07	15
Bepridil ^a	0.96	-ve inotrope	0.7	22
Chlorpromazine ^a	1.04	-ve inotrope	1.02	28
Cisapride ^a	0.26	-ve inotrope	0.02	8
Clarithromycin ^a	120	-ve inotrope	16	13
Clozapine ^a	2.13	-ve inotrope	1.5	21
D, L-Sotalol ^a	450	No effect	>450	>30
Disopyramide ^a	21	-ve inotrope	9.3	13
Dofetilide ^a	0.2	No effect	>0.2	>100
Domperidone ^a	2	-ve inotrope	0.2	10
Droperidol ^a	0.48	-ve inotrope	0.18	11
Erythromycin	5.1	No effect	>5.1	>30
Flecainide	22.6	-ve inotrope	1.1	2
Ibutilide ^a	3	-ve inotrope	2	20
Moxifloxacin	329	No effect	>329	>30
Ondansetron ^a	11.2	-ve inotrope	14	34
Procainamide	1625	-ve inotrope	2215	38
Quinidine ^a	100	-ve inotrope	3.6	1
Sematilide	133	No effect	>133	>30
Terodiline	4.35	-ve inotrope	0.7	5
Vandetanib ^a	9	-ve inotrope	2.7	9
Diltiazem ^a	3.84	-ve inotrope	1	8
Diphenhydramine	1.02	-ve inotrope	0.6	17
Loratadine ^a	0.0135	-ve inotrope	0.0175	35
Mexiletine ^a	75	-ve inotrope	0.9	0.4
Mibefradil	0.36	-ve inotrope	0.18	13
Nifedipine ^a	0.23	-ve inotrope	0.04	5
Nitrendipine ^a	0.091	-ve inotrope	0.06	18
Ranolazine ^a	200	-ve inotrope	17	9
Tamoxifen ^a	0.663	-ve inotrope	0.99	36
Verapamil ^a	10	-ve inotrope	0.04	2

IC₅₀: Concentration inducing 50% decrease in sarcomere shortening; Hill equation using SigmaPlot v13 was fitted to sarcomere shortening concentration-effect curves, assuming drugs would eventually cause complete inhibition of the contractility when they decreased sarcomere shortening by ≥25%. ^a: CiPA-selected drug; fETPC, free effective therapeutic plasma concentration.

Nguyen et al., 2017 FiP

Inhibition of Kinase Activity to Control Tumor Growth


Can Lead to Cardiotoxicity




Validating Clinical Relevance of Cancer Agents

- Validated 9 clinical well characterized controls :
 - 1) 4 toxic TKIs (Sorafenib, Vandetanib, AZD7762, Imatinib)
 - 2) 4 non-toxic TKIs (Erlotinib, Dasatinib, Afatinib, Gefitinib)
 - 3) One toxic anthracycline (Doxorubicin)
 - 4) Each drug was tested at multiples of the Cmax
 - 5) Each concentration was perfused for 5 mins

Sorafenib Induces Functional Cardiotoxicity




Afatinib Induces No Functional or Structural Cardiotoxicity



*0.3 μ M = 3-fold C_{max}
1 μ M = 10-fold C_{max}
3 μ M = 30-fold C_{max}*

Tyrosine Kinase Inhibitors Affect Adult Human Primary Cardiomyocyte Contractility



TKI	Clinical contractility risk	Human cardiomyocyte contractility	C _{max} (μM)	IC ₅₀ (μM)	Ratio (IC ₅₀ /C _{max})
Sorafenib			3.4	1.2	0.35
Vandetanib			1.8	4.6	2.55
AZD7762			0.12	0.8	6
Imatinib			5	44	8
Erlotinib			2.5	>10*	>4
Dasatinib			0.01	>0.3	>30
Afatinib			0.1	>3	>30
Gefitinib			0.1	>3	>30


Similar human cardiac tissue data recently published by Schneider C et al., 2018 Nature Scientific Reports

*: Limit of solubility

AZD7762 Heart Data Case study


AZD7762 Animal Heart Data

No Preclinical to Clinical Translation







AZD7762 Human Heart Data


Preclinical to Clinical Translation



Low Inter- and Intra-Heart Variability



Safety Index and Expression of TKI Pro-arrhythmia Risk




TKI	Clinical contractility risk	Human cardiomyocyte pro-A	C _{max} (μM)	Pro-A conc. (μM)	Ratio (Pro-A conc./C _{max})
Vandetanib	Risk	Risk	1.8	3	1.6
AZD7762	No risk	Risk	0.12	0.3	2.5
Imatinib	Risk	Risk	5	6	1.2
Sorafenib	No risk	No risk	3.4	>10	>3
Erlotinib	No risk	No risk	2.5	>10*	>4
Dasatinib	No risk	No risk	0.01	>0.3	>30
Afatinib	No risk	No risk	0.1	>3	>30
Gefitinib	No risk	No risk	0.1	>3	>30

*: Limit of solubility; Pro-A: Pro-arrhythmia



■ Risk
■ No risk

Doxorubicin, Anthracycline Agent, Affects Adult Human Primary Cardiomyocyte Contractility



Safety Index and Expression of Doxorubicin

Pro-arrhythmia Risk




Positive Inotropy Assessment - Validation Set

Targets 12 Different Mechanisms of Action

Mechanism of Action	Drug
Na ⁺ /K ⁺ pump inhibition	Digoxin
Na ⁺ /K ⁺ pump inhibition	Ouabain
Na ⁺ /Ca ²⁺ exchanger inhibition	SEA-0400
Myosin activation	Omecamtiv Mecarbil
Myosin activation	EMD-57003
Ca ²⁺ sensitization	Levosimendan
Non-selective b-adrenoceptor activation	Isoproterenol
Non-selective b-adrenoceptor activation	Epinephrine
b1-adrenoceptor activation	Dobutamine
PDE3 inhibition	Milrinone
PDE inhibition	IBMX
Ca ²⁺ channel activation	Bay-K 8644
Adenylyl cyclase activation	Forskolin
Adenylyl cyclase activation	NKH-477
SERCA activation	N106
RyR activation	Caffeine










Isoproterenol Induces Positive Inotropic Effect




Baseline level

Iso level


Identification of Positive Inotropic Effects and Determination of Exposure Responses



Differential Effects of Positive Inotropes




Cluster Analysis Is Used To Fingerprint Compounds With Inotropic Effects




Heatmap data generated from 4th concentration data. Red and green colors indicate decrease and increase of >25% and 10% change, respectively. Black colors indicate no effect (<-25% < % change < 10%). Numbers in boxes indicate means % change relative to vehicle.

Contraction Failure and Pause-Dependent Arrhythmia

Compounds With Na⁺ Channel Liability




Segregation of Ca^{2+} -Dependent Mechanisms



2D PCA generated from top test concentration data. Blue and red colors indicate increase and decrease in Ca^{2+} , respectively. Ellipses show confidence intervals of 0.75.

Adult Human Primary Cardiomyocyte Model

Integrated Into the CiPA Paradigm




This new potential approach was presented to the ICH S7B & E14 Discussion Group Meeting (Charlotte, NC, November 13, 2018).

Higher Throughput Technology for Adult Human Primary Cardiomyocytes


- Simultaneous sarcomere shortening measurement from multiple myocytes
- Automated acquisition & analysis
- Same drug exposure time for all cells (minimum of 5 cells)
- Screening capacity: 10 to 40 compounds / per day
- Simultaneous contractility & pro-arrhythmia assessment

MyoBLAZER™ Proprietary Technology

A





B



Adult Human Primary Cardiomyocyte Model

Early Primary Cardiotoxicity Detection Tool



Adult Human Primary Cardiomyocyte Model

Early Primary Cardiotoxicity Detection Tool

- Permit rigorous and integrated human cardiac drug discovery at the preclinical stages
- Differentiate safe from cardiotoxic drugs
- Can enable mechanistic assessment
- Predictive of clinical outcomes

AnaBios Ex-Vivo Adult Human Primary Platforms

Future Assays and Products

- Adult primary cardiomyocytes
 - ✓ Cell preservation & storage
 - ✓ Cell shipment & distribution
 - ✓ Distribution of MyoBLAZER™
- Human lung slices
 - ✓ Normal
 - ✓ Diseased fibrotic
- Diseased human hearts
 - ✓ LVH
 - ✓ CAD / MI
 - ✓ HF
 - ✓ Afib
 - ✓ Diabetic
- Human liver slices
 - ✓ Normal
 - ✓ NASH
 - ✓ Steatosis

AnaBios - SPS 2019 Meeting

Booth 503 & Posters

- Poster 023: Predicting Cardiotoxicity of Cancer Tyrosine Kinase Inhibitors with Adult Human Primary Cardiomyocytes
 1. Tuesday, September 24, 15:15-15:45
 2. Wednesday, September 25, 13:00-14:00
- Poster 098: Late Sustained Sodium Current in Adult Human Primary Cardiomyocytes
 1. Tuesday, September 24, 13:00-14:00
 2. Wednesday, September 25, 15:15-15:45

Thank You!