


Introduction

The late sustained sodium current ($I_{Na,L}$), a depolarizing current that persists throughout the action potential (AP) plateau, contributes to the AP duration and maintains the intracellular homeostasis of Na^+ . Increase or inhibition of $I_{Na,L}$ is often associated with arrhythmogenicity or mitigation of pro-arrhythmia risk, respectively. Indeed, an increased $I_{Na,L}$ has been associated with the long QT syndrome type 3 (LQT 3). Since $I_{Na,L}$ was one of the selected channels for the CiPA initiative and drugs that block the hERG channel and also inhibit $I_{Na,L}$ are not associated with pro-arrhythmia in humans, identifying the effect of compounds on $I_{Na,L}$ in human cardiomyocytes during preclinical development can aid in the determination of pro-arrhythmia risk for new drugs.

Methods


In order to better understand the properties of this current, we performed voltage-clamp recording of $I_{Na,L}$, at physiological temperature, in whole-cell patch-clamp experiments using adult human primary cardiomyocytes isolated from ethically consented donor hearts. Test articles were applied until steady state effect was achieved.

$I_{Na,L}$ CiPA voltage-clamp recording protocol




$I_{Na,L}$ protocol was elicited every 10 s from a holding potential of -80 mV. $I_{Na,L}$ current was measured during the step-ramp, which goes from 40 mV to -95 mV over a duration of 100 ms.

Ranolazine inhibits and ATX II stimulates $I_{Na,L}$



A) Myocyte was first exposed to control solution. Applications of Ranolazine and ATX II are indicated by the lines.
B) Typical $I_{Na,L}$ traces. The dotted line indicates the zero-current level.
C) $I_{Na,L}$ current density-voltage relationship ($n=3$ cells).

Ranolazine inhibits the ATX II stimulation of $I_{Na,L}$



A) Myocyte was first exposed to control solution. Perfusion with ATX II alone or in combination with Ranolazine are indicated by the lines.
B) Typical $I_{Na,L}$ traces. The dotted line indicates the zero-current level.
C) $I_{Na,L}$ current density-voltage relationship ($n=6$ cells).

Summary

1. Human cardiomyocytes express functional $I_{Na,L}$ and can differentiate $I_{Na,L}$ inhibitors from facilitators.
2. Action potential experiments are currently underway to assess $I_{Na,L}$ role in modulating pro-arrhythmia.
3. Human cardiomyocytes could potentially provide a useful strategy for the early assessment of the ability of new multichannel blocking drugs with $I_{Na,L}$ affinity to prevent the occurrence of pro-arrhythmia.