
B370

Introduction

The late sustained sodium current (INa,L), a depolarizing current that persists throughout the action potential (AP) plateau, contributes to the AP duration and maintains the intracellular homeostasis of Na⁺. Increase or inhibition of INa,L is often associated with arrhythmogenicity or mitigation of pro-arrhythmia risk, respectively. Since INa,L was one of the selected channels for the CiPA (Comprehensive In Vitro Pro-arrhythmia Assay) initiative and drugs that block the hERG channel and inhibit INa,L are not associated with pro-arrhythmia in humans, identifying the effect of compounds on INa,L in human cardiomyocytes during preclinical development can aid in the determination of pro-arrhythmia risk for novel drugs.


Methods

In order to better understand the properties of INa,L, we performed voltage-clamp and current-clamp recordings at physiological temperature (37°C), in whole-cell patchclamp experiments using adult human primary ventricular cardiomyocytes isolated from ethically consented donor hearts. INa,L enhancer (ATX-II) and inhibitors (Ranolazine and GS967) were applied until a steadystate effect was achieved.



INa,L Voltage-Clamp Recording Protocol

AnaBios Corp., San Diego, CA 92109, USA



simultaneous perfusion group. Solid black line: 0 mV level.